

## PROTOCOL: Illumina Barcoded Paired-End Whole Genome Shotgun Library Preparation

This protocol provides instructions for preparing barcoded paired-end whole genome shotgun (WGS) libraries for sequencing by Illumina platforms (GAII, HiSeq and MiSeq). It involves using the Covaris S2 system for shearing DNA samples, using the NEBNext End Repair, A-Tailing, and Ligation Modules for DNA modification, as well as using the 2X Phusion High Fidelity PCR Master Mix for ligation-mediated PCR (LM-PCR). Quality control standard or criteria for each major procedural step is also addressed here.

The whole process described in this protocol includes the following steps:

### WGS Library preparation

DNA fragmentation

DNA end repair

3'-end adenylation

Ligation with Illumina multiplexing PE adaptor

LM-PCR, purification and QC

### Reagents

Item
NEBNext End-Repair Module, Cat# E6050L, NEB
NEBNext dA-Tailing Module, Cat# E6053L, NEB
NEBNext Ligation Module, Cat# E-6056L, NEB
<b>Illumina Multiplexing PE Adaptor Oligos</b> Note: All these primer oligos need to be HPLC purified.  Illumina multiplexing PE adaptor oligo A 5'-ACACTCTTTCCCTACACGACGCTCTTCCGATCT  Illumina multiplexing PE adaptor oligo B 5'-phospho-GATCGGAAGAGCACACGTCT
<b>LM-PCR primer IMUX-P1.0 (working concentration 50uM)</b> 5'-AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATC*T  HPLC purified; *addition of phosphothioate bond.

**LM-PCR barcoded primers (working concentration 50uM)**

Note: All these primer oligos need to be HPLC purified. Addition of **phosphothioate bond** before addition of the last "T" as a preventative measure against degradation in storage.

**6-bp barcode set (the barcodes were designed by Illumina Inc.)**

IBC1 CAAGCAGAAGACGGCATAACGAGATCGTGATGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IBC2 CAAGCAGAAGACGGCATAACGAGATACATCGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IBC3 CAAGCAGAAGACGGCATAACGAGATGCCTAAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IBC4 CAAGCAGAAGACGGCATAACGAGATTGGTCAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IBC5 CAAGCAGAAGACGGCATAACGAGATCACTGTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IBC6 CAAGCAGAAGACGGCATAACGAGATATTGGCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IBC7 CAAGCAGAAGACGGCATAACGAGATGATCTGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IBC8 CAAGCAGAAGACGGCATAACGAGATTCAAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IBC9 CAAGCAGAAGACGGCATAACGAGATCTGATCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IBC10 CAAGCAGAAGACGGCATAACGAGATAAGCTAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IBC11 CAAGCAGAAGACGGCATAACGAGATGTAGCCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IBC12 CAAGCAGAAGACGGCATAACGAGATTACAAGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT

**9-bp barcode set (the barcodes were designed by Adey et. al , Genome Biol. 2010; 11: R119)**

IDMB1 CAAGCAGAAGACGGCATAACGAGATTACGAAGTCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB2 CAAGCAGAAGACGGCATAACGAGATGACGAGATTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB3 CAAGCAGAAGACGGCATAACGAGATACCGTAAGAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB4 CAAGCAGAAGACGGCATAACGAGATTAGTGGAAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB5 CAAGCAGAAGACGGCATAACGAGATCATTACCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB6 CAAGCAGAAGACGGCATAACGAGATTCGTTGAAGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB7 CAAGCAGAAGACGGCATAACGAGATTAGTACGCTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB8 CAAGCAGAAGACGGCATAACGAGATCTCAGATCAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB9 CAAGCAGAAGACGGCATAACGAGATTTACCGTAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB10 CAAGCAGAAGACGGCATAACGAGATGTCATGCATGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB11 CAAGCAGAAGACGGCATAACGAGATAGGACAGTTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB12 CAAGCAGAAGACGGCATAACGAGATATGGTGTCTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB13 CAAGCAGAAGACGGCATAACGAGATGGATGTTCTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB14 CAAGCAGAAGACGGCATAACGAGATCTTATCCAGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB15 CAAGCAGAAGACGGCATAACGAGATGTAAGTCACGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB16 CAAGCAGAAGACGGCATAACGAGATTTCAAGTGAGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB17 CAAGCAGAAGACGGCATAACGAGATCTCGTAATGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB18 CAAGCAGAAGACGGCATAACGAGATCATGTCTCAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB19 CAAGCAGAAGACGGCATAACGAGATAATCGTGGAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB20 CAAGCAGAAGACGGCATAACGAGATGTATCAGTCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB21 CAAGCAGAAGACGGCATAACGAGATAGCAGATGTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB22 CAAGCAGAAGACGGCATAACGAGATTCCTAACGTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB23 CAAGCAGAAGACGGCATAACGAGATAACAGTCCAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB24 CAAGCAGAAGACGGCATAACGAGATCCTTGAGAAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT  
 IDMB25 CAAGCAGAAGACGGCATAACGAGATTTAAGCCTGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT

IDMB26	CAAGCAGAAGACGGCATAACGAGATTTAGACCACGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB27	CAAGCAGAAGACGGCATAACGAGATTGTCTAGTGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB28	CAAGCAGAAGACGGCATAACGAGATTAGATCGAGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB29	CAAGCAGAAGACGGCATAACGAGATTGAATGCCAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB30	CAAGCAGAAGACGGCATAACGAGATGTGCAATGTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB31	CAAGCAGAAGACGGCATAACGAGATAGTGGCATAAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB32	CAAGCAGAAGACGGCATAACGAGATATGATCGGTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB33	CAAGCAGAAGACGGCATAACGAGATAGTCTACCTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB34	CAAGCAGAAGACGGCATAACGAGATGATCAACTGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB35	CAAGCAGAAGACGGCATAACGAGATATCGGTAGTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB36	CAAGCAGAAGACGGCATAACGAGATCGTATGATGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB37	CAAGCAGAAGACGGCATAACGAGATTTACTGACGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB38	CAAGCAGAAGACGGCATAACGAGATCTGTCGTAAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB39	CAAGCAGAAGACGGCATAACGAGATCAACTGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB40	CAAGCAGAAGACGGCATAACGAGATATCGATCTCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB41	CAAGCAGAAGACGGCATAACGAGATGCAACTATGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB42	CAAGCAGAAGACGGCATAACGAGATGATGACTTCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB43	CAAGCAGAAGACGGCATAACGAGATGACGTTACAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB44	CAAGCAGAAGACGGCATAACGAGATCATCTGCTAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB45	CAAGCAGAAGACGGCATAACGAGATATTAGTCGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB46	CAAGCAGAAGACGGCATAACGAGATTAGCGTACTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB47	CAAGCAGAAGACGGCATAACGAGATCCAAGCAATGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB48	CAAGCAGAAGACGGCATAACGAGATCCGTAATTGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB49	CAAGCAGAAGACGGCATAACGAGATAGAATTGCCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB50	CAAGCAGAAGACGGCATAACGAGATACCTGTAACGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB51	CAAGCAGAAGACGGCATAACGAGATCATCAGTGTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB52	CAAGCAGAAGACGGCATAACGAGATGAATCCTCAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB53	CAAGCAGAAGACGGCATAACGAGATGCTGTATACGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB54	CAAGCAGAAGACGGCATAACGAGATGAAGGCTATGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB55	CAAGCAGAAGACGGCATAACGAGATGGAATCGATGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB56	CAAGCAGAAGACGGCATAACGAGATGCTTATGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB57	CAAGCAGAAGACGGCATAACGAGATTGACGCATTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB58	CAAGCAGAAGACGGCATAACGAGATCACGATTCTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB59	CAAGCAGAAGACGGCATAACGAGATTATTGCCTCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB60	CAAGCAGAAGACGGCATAACGAGATAAGTCAGAGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB61	CAAGCAGAAGACGGCATAACGAGATATAGCTGAGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB62	CAAGCAGAAGACGGCATAACGAGATTGCTCACAAAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB63	CAAGCAGAAGACGGCATAACGAGATGTCTTCTGAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB64	CAAGCAGAAGACGGCATAACGAGATTTGCCGATTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB65	CAAGCAGAAGACGGCATAACGAGATCTCGAATACGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB66	CAAGCAGAAGACGGCATAACGAGATTGGCTTCTAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB67	CAAGCAGAAGACGGCATAACGAGATAAGGCCATTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB68	CAAGCAGAAGACGGCATAACGAGATAAGTTGACCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT

IDMB69	CAAGCAGAAGACGGCATAACGAGATCTGAACTGAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB70	CAAGCAGAAGACGGCATAACGAGATCTAGGTGTAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB71	CAAGCAGAAGACGGCATAACGAGATCCATCTTAGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB72	CAAGCAGAAGACGGCATAACGAGATCTACGACATGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB73	CAAGCAGAAGACGGCATAACGAGATTCCAACATGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB74	CAAGCAGAAGACGGCATAACGAGATGCTATCATCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB75	CAAGCAGAAGACGGCATAACGAGATACAGCTTCAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB76	CAAGCAGAAGACGGCATAACGAGATAGTCATTGCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB77	CAAGCAGAAGACGGCATAACGAGATAGATCTCGAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB78	CAAGCAGAAGACGGCATAACGAGATATGCTCTTGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB79	CAAGCAGAAGACGGCATAACGAGATTTAGTGCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB80	CAAGCAGAAGACGGCATAACGAGATTCTAGTTCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB81	CAAGCAGAAGACGGCATAACGAGATGGTGCATTAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB82	CAAGCAGAAGACGGCATAACGAGATACTGAGGATGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB83	CAAGCAGAAGACGGCATAACGAGATTAGCAGTCAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB84	CAAGCAGAAGACGGCATAACGAGATTCAGTCTCGAAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB85	CAAGCAGAAGACGGCATAACGAGATAACCAATCAGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB86	CAAGCAGAAGACGGCATAACGAGATGATATGGACGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB87	CAAGCAGAAGACGGCATAACGAGATTGAGAGATCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB88	CAAGCAGAAGACGGCATAACGAGATTGCCATTAGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB89	CAAGCAGAAGACGGCATAACGAGATACTAACGCAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB90	CAAGCAGAAGACGGCATAACGAGATATGTAGCACGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB91	CAAGCAGAAGACGGCATAACGAGATGGTCGATATGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB92	CAAGCAGAAGACGGCATAACGAGATGCGAGTTATGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB93	CAAGCAGAAGACGGCATAACGAGATGACTGAGTAGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB94	CAAGCAGAAGACGGCATAACGAGATAGATACTCCGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB95	CAAGCAGAAGACGGCATAACGAGATGCTAGAGTTGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
IDMB96	CAAGCAGAAGACGGCATAACGAGATAATGTAGCGGTGACTGGAGTTCAGACGTGTGCTCTTCCGATCT
Agencourt AMPure XP beads(60ml), Cat# A63882, Beckman Coulter	
2X Phusion High Fidelity PCR Master Mix with HF buffer, Cat# M0531, NEB	

## Consumables

Item	Source	Catalog #
Agilent Chip 7500	Agilent	5067-1506
2.2% FLASH DNA gel cassette	Lonza/VWR	57023
Flash Gel Loading Dye (5X)	Lonza/VWR	50462
Flash Gel Marker 100 – 4kb	Lonza/VWR	50473
EtOH, 200 proof	Aaper	Stock room
Clear 12-strip 0.2ml PCR tubes w/ caps	Phenix	MPC-665
Caps Only 12-Strip Clear for MPC-655	Phenix	MPC-660

1.7 ml microtubes (MCT-175-B)	VWR	10011-720
10 ul Pipet tips	VWR	47449-868
200 ul pipet tips	VWR	22234-016
1000 ul pipet tips	VWR	16466-008
Shearing microtube	Covaris	520045
Nuclease-free Water	Qiagen	129114
Flashgel Starter Kit	VWR	95015-612

### Equipments and Devices

Item	Source	Catalog #
Covaris S2 Acoustic Shearer	Covaris	600028
P1000 pipette	Rainin	L-1000
P200 pipette	Rainin	L-200
P10 pipette	Rainin	L-10
Gel Imager PC 2000	BioRad	1708195
Thermomixer	Eppendorf	0022670000
Agilent Bioanalyzer 2100	Agilent	G2940CA
Microcentrifuge 5424	Eppendorf	5424000410
AB GeneAmp Thermocycler	Applied Biophysics	N8050200
Thermocycler Sample Block Module	Applied Biophysics	4314443
DynaMag2 particle concentration rack	Invitrogen	123-21D

#### **Before you start, make sure that you read carefully the following safety information:**

- Wear lab coat, gloves, and protective goggles.
- Protect bench working area with absorbent bottom plasticized pad.
- Discard all materials in Biohazard bag. Discard all liquid in a Biohazard labeled plastic bottle that should be tightly closed and disposed of in the Biohazard bag when full.

#### **References:**

- Illumina Multiplexing\_SamplePrep\_Guide\_1005361\_D.


## **WGS Library Preparation**

### **1. DNA Fragmentation**

- a) Fill a Covaris™ water bath to the level of 12 in the Covaris S2 device and degas for 30-45 minutes prior to shearing your gDNA in the S2. When a microtube is placed in the tube holder, the water level should be at the base of the cap and the glass portion of the tube should be completely submerged. The water bath temperature should be between 6 °C and 8 °C.

- b) Dilute 500ng of DNA (based on pico-green or Qubit) into 100  $\mu$ l of low-TE water in a microcentrifuge tube.
- c) Place a Covaris™ microTube (6X16mm) into the loading station. Keep the cap on the tube and use a tapered pipette tip to slowly transfer the 100  $\mu$ L of DNA sample through the pre-split septa. Be careful not to introduce a bubble into the bottom of the tube. Bubbles will interfere with the acoustic field.
- d) Shear the DNA using the following Covaris™ S2 System shearing conditions.
  1. Number of Cycles: **2 (40 seconds per cycle)** (if sample quality is poor, start with 1 cycle and check on a gel).
  2. Batch Temperature: **6-8 °C**
  3. Bath Temperature Limit: **<15 °C**
  4. Mode: **Frequency sweeping**
  5. Water Quality Testing Function: **Off**
  6. Duty cycle: **10%**
  7. Intensity: **4**
  8. Cycles/burst: **200**
  9. Time: **80 s**
- e) Place the Covaris™ microtube into the loading station. While keeping the snap-cap on, insert a pipette tip through the pre-split septa, slowly remove the sheared DNA and transfer it into a new 1.7 ml tube.
- f) Load 1  $\mu$ l of each sample and 2  $\mu$ l of 100 bp – 4000 bp marker on a 2.2% FlashGel DNA Cassette (Lonza, Cat#: 57023). Gel running condition set as the following: Voltage 270, running time: 4 minutes. Gel imager (Bio-Rad) setting should be uniformed to the condition as described here. Condition of Bio-Rad gel imager: Filter 1 for Eth Bromide, auto expose of 0.6 s. Image transform: Higher: 3000, low: 200, Gamma: 1.00, Min: 79, Max: 4095.

**Note:** The average size should be approximately 300 bp and the majority of fragment sizes should be smaller than 500 bp (see Figure 1 gel picture for example). If majority size of sheared fragments is larger than 800 bp, re-shear the sample with 20 s increment each time for up to 2 additional cycles. The negative control should give no product.



**Figure 1 Post-Covaris DNA size distribution on 2.2% Flash Gel. Marker lane: FlashGel DNA marker (100 bp – 4000 bp, Cat. No. 50473, Lonza).**

- g) Purify by adding 1.8X Agencourt AMPure XP beads and binding by rotating on a thermomixer for about 5 minutes at RT before placing into the magnetic particle concentrator (MPC). Remove the supernatant (and save). Wash 3 times with freshly prepared 70% ethanol. Keep the sample tube in MPC during the wash. Add ethanol slowly to avoid disturbing the bead pellet. After wash, briefly dry the beads in the air. Elute DNA in 52  $\mu$ l of nuclear-free H<sub>2</sub>O.

## 2. DNA End Repair

- a) Combine and mix the following components in a 1.7 ml tube:

Component	Volume ( $\mu$ L)
Sheared DNA	50.0
End-repair 10x buffer*	9.0
End-repair enzyme mix*	5.0
Nuclease-free water	26.0
Total	90.0

\*From NEBNext End-Repair Module (Cat. No. E6050L).

- b) Incubate the mixture at 20°C for 30 minutes at bench top cooler.
- c) Purify by adding 1.8X Agencourt AMPure XP beads and binding by rotating on a thermomixer for about 5 minutes at RT before placing into the magnetic particle concentrator (MPC). Remove the supernatant (and save). Wash 3 times with freshly prepared 70% ethanol. Keep the sample tube in MPC during the wash.

Add ethanol slowly to avoid disturbing the bead pellet. After wash, briefly dry the beads in the air. Elute DNA in 42  $\mu$ l of nuclear-free H<sub>2</sub>O.

### 3. 3'-end Adenylation

a) Combine and mix the following components in a 1.7 ml tube:

Component	Volume ( $\mu$ L)
End-repaired DNA	40.0
NEBNext <sup>TM</sup> dA-Tailing Reaction Buffer (10X)*	6.0
Klenow Fragment (3'-5' exo <sup>-</sup> )*	3.0
Nuclease-free water	11.0
Total	60.0

\*From NEBNext dA-Tailing Module (Cat. No. E6053L).

b) Incubate the mixture at 37°C for 30 min in a calibrated thermomixer.

4) Purify by adding 1.8X Agencourt AMPure XP beads and binding by rotating on a thermomixer for about 5 minutes at RT before placing into the magnetic particle concentrator (MPC). Remove the supernatant (and save). Wash 3 times with freshly prepared 70% ethanol. Keep the sample tube in MPC during the wash. Add ethanol slowly to avoid disturbing the bead pellet. After wash, briefly dry the beads in the air. Elute DNA in 42  $\mu$ l of nuclear-free H<sub>2</sub>O.

### 4. Ligate Adaptors to the DNA

Oligonucleotide annealing:

**Illumina multiplexing PE adapter oligo A and Illumina multiplexing PE adapter oligo B should be ordered as HPLC purified.**

Mix the paired oligonucleotides at final concentration of 300  $\mu$ M in 1x Ligase buffer (from NEBNext Ligation Module, Cat# E6056-L) and run with the following annealing program on thermocycler: 95°C for 5 min, 80°C for 3 min, 70°C for 3 min, 60°C for 3 min, 50°C for 3 min, 40°C for 3 min, 30°C for 3 min, 20°C for 3 min and 4°C hold. Dilute 20x to obtain 15  $\mu$ M working concentration. Aliquot annealed adaptor into eppendorf tubes and store them in -20°C freezer.

a) Combine and mix the components below

Component	Volume ( $\mu$ L)
-----------	-------------------



Nuclease-free water	22.0
Illumina Multiplexing PE Adaptors (15 $\mu$ M concentration)	5.0
Quick Ligase 5X buffer*	18.0
A-Tailed DNA	40.0
Quick Ligase Enzyme*	5.0
Total	90.0

- From NEB (Cat. No. E-6056L).

- b) Incubate at room temperature for 30 minutes.
- c) Purify by adding 1.0X Agencourt AMPure XP beads and binding by rotating on a thermomixer for about 5 minutes at RT before placing into the magnetic particle concentrator (MPC). Remove the supernatant (and save). Wash 3 times with freshly prepared 70% ethanol. Keep the sample tube in MPC during the wash. Add ethanol slowly to avoid disturbing the bead pellet. After wash, briefly dry the beads in the air. Elute DNA in 58  $\mu$ l of nuclease-free H<sub>2</sub>O.

**Note:** (Optional) Assess yield and quality by running Agilent Bioanalyzer 2100 DNA 7500 Chip. Make sure there is no excess free adaptor's peak in the eluted DNA sample.

## 5. LM-PCR Amplification


- a) Prepare the following PCR mix: using only A or B column instructed as follow for Invitrogen enzyme mix or NEB phusion enzyme mix accordingly

Component	Volume ( $\mu$ L)
DNA ligation product	56.0
2X Phusion High Fidelity PCR Master Mix with HF buffer	60.0
LM-PCR primer IMUX-P1.0 (50 $\mu$ M)	2.0
LM-PCR barcoded primer IBCn (n=1-12) or IDMBn (n=1-96) (50 $\mu$ M)	2.0
Total	120.0

- b) Vortex gently and transfer mix into a 0.2 ml PCR strip tubes. Seal the wells with caps.
- c) Place them in the ABI GeneAmp PCR System 9700/Veriti for amplification enrichment. PCR condition:

(1) 30" @ 98°C

- (2) 10" @ 98°C
  - (3) 30" @ 60°C
  - (4) 30" @ 72°C
  - (5) Repeat step (2) to (4) 6 times for total 6 cycles
  - (6) 5' @ 72°C
  - (7) HOLD @ 4°C
- d) Load 1 µl of sample into one lane of a 2.2% FlashGel to check amplification.
- e) Purify by adding 1.0X Agencourt AMPure XP beads and binding by rotating on a thermomixer for about 5 minutes at RT before placing into the magnetic particle concentrator (MPC). Remove the supernatant (and save). Wash 3 times with freshly prepared 70% ethanol. Keep the sample tube in MPC during the wash. Add ethanol slowly to avoid disturbing the bead pellet. After wash, briefly dry the beads in the air. Elute DNA in 40 µl of nuclear-free H<sub>2</sub>O.
- f) Run an Agilent Bioanalyzer 2100 DNA Chip 7500 to check size distribution and quantity of PCR product. Make sure there is no excess primer left in the eluted sample (See Figure 2 for example). The yield should be at least >1µg.



**Figure 2 Example of Illumina barcoded paired end WGS library on Agilent Bioanalyzer 2100 DNA 7500 Chip.**

**(end)**